cefar[®] peristim^{pro}

	USERS GUIDE	3
	MANUAL DEL USUARIO	25
	GEBRAUCHSANWEISUNG	47
	MODE D'EMPLOI	69
	ISTRUZIONI PER L'USO	91
	GEbruikersHANDLEIDING	113
	BRUKSANVISNING	135

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OPERATING INSTRUCTIONS – ENGLISH

BEFORE USING THE STIMULATOR

1. INTRODUCTION	4
2. MEDICAL BACKGROUND	4
3. PRECAUTIONARY MEASURES	6

USER INFORMATION

4. OVERVIEW	10
CONTROL BUTTONS	10
DISPLAY SYMBOLS	11
5. OPERATION	12
TIMER	14
PROGRAM PAUSE	14
6. REPLACEMENT OF BATTERIES	15
7. PROGRAMS	16
PRESET PROGRAMS	16
CUSTOM PROGRAMS (P10-P12)	18
PROGRAM LOCK	20
COMPLIANCE	20

ADDITIONAL INFORMATION

8. ACCESSORIES	21
9. CARE INSTRUCTIONS	22
10. TROUBLESHOOTING	23
11. TECHNICAL DATA	24
KEY TO THE SYMBOLS	25
12. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES	26
13. REFERENCES	30

The Cefar stimulators are widely used by public and private health care professionals around the world.

Electrical nerve and muscle stimulation is effective, has no side effects and it is economical. Through clinical research, areas of application for TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscle Stimulation) are rapidly expanding. Cefar is working actively to further develop the method towards a natural treatment alternative for both health care professionals and consumers.

More information about TENS, EMS and our products can be found on our web site.

INTENDED USE:

CEFAR PERISTIM PRO (1527) is a dual channel incontinence stimulator for treatment and rehabilitation of incontinence. CEFAR PERISTIM PRO features seven preset incontinence program for internal use with a vaginal and/or anal probe, two preset TENS program for external use with two surface electrodes. The CEFAR PERISTIM PRO provides programming for three custom programs. The channels are simultaneous, which means that a selected program applies for both channels.

INTENDED USER:

The user of CEFAR PERISTIM PRO can be a Health Care Professional or a patient. The device should be used indoors and may be used in a healthcare facility setting or in a home environment.

2. MEDICAL BACKGROUND

INCONTINENCE

Urinary incontinence, involuntary loss of urine from the bladder, is a problem for many people. There are two main types of urinary incontinence; stress incontinence and urge incontinence. Faecal incontinence, involuntary passage of faeces, is not often discussed, but still a common problem. Electrical stimulation through a vaginal/anal probe, or in some cases surface electrodes, is a well-tolerated treatment for urge, stress, mixed and faecal incontinence and has shown positive results in improving bladder and bowel control.

Stress incontinence

Stress incontinence is urine leakage caused by increased abdominal pressure on the bladder, such as coughing, sneezing, laughing, exercising or lifting something heavy. Stress incontinence is the most common type of incontinence and is primarily affecting women. It usually occurs when the perineal and pelvic floor muscles are weakened, for example by pregnancy, childbirth or menopause.

Urge incontinence

Urge incontinence means a sudden, strong urge to urinate followed by an immediate bladder contraction, resulting in an involuntary leakage of urine. Both men and women can be affected, particularly the elderly. One reason for this condition can be a disruption in the part of the nervous system that controls the bladder.

Mixed incontinence

Mixed incontinence is a combination of stress incontinence and urge incontinence.

Faecal incontinence

Faecal incontinence, also called anal or bowel incontinence, is the impaired ability to control passage of gas or stool. There are many possible causes of faecal incontinence, the most common is injury to the anal sphincter (ring-like muscle), for instance during childbirth or surgery, or damage to the nerves that control the anal sphincters. The condition usually becomes worse as people age.

INCONTINENCE TREATMENT

Electrical stimulation via the pelvic nerves is a recognized treatment alternative for urinary incontinence. It is also proposed as a treatment method for faecal incontinence due to pelvic floor dysfunction or a poorly functioning anal sphincter.

When treating stress incontinence, the aim of the electrical stimulation is to mirror a voluntary muscle contraction and improve the function of the pelvic floor muscles. For urge incontinence, the aim is to inhibit involuntary bladder contractions by stimulating the nerves in the pelvic floor. When treating mixed incontinence, a stimulation appropriate for both urge incontinence and stress incontinence is used. For faecal incontinence, the aim is to improve bowel control by strengthening and toning pelvic floor muscles.

- ## WARNINGS

- People with implanted electronic equipment, such as pacemakers and intracardiac defibrillators, must not be treated with CEFAR PERISTIM PRO.
- Pregnant women must not be treated with the urology programs.
- People with extra-urethral incontinence (fistula, ectopic ureter) must not be treated with the urology programs.
- People with overflow incontinence due to outflow obstacle must not be treated with the urology programs.
- People with serious retention of urine in the upper urinary tract must not be treated with the urology programs.
- People with complete peripheral denervation of the pelvic floor must not be treated with the urology programs.
- Stimulation should not take place while the user is connected to high-frequency surgical equipment. It may cause burn injuries on the tissue, as well as problems with the stimulator.
- Do not use the stimulator in the vicinity of shortwave or microwave therapy equipment, since this may affect the output power of the stimulator.
- Keep the stimulator out of reach of children.

- Patients with total/subtotal prolapsed uterus/vagina should be stimulated with greatest caution.
- Patients with urinary tract infections must be treated and clear of infection before starting therapy with the urology programs. Consult your doctor.
- If tissue irritation should occur, treatment should be temporarily discontinued. If problems continue, contact your health care provider. Hypersensitivity can occur in isolated cases. The problem usually disappears when probes or gel are changed to another type.
- Observe caution when using electrotherapy at the same time as the patient is connected to monitoring equipment with body worn electrodes. The stimulation might interfere with the signals to the monitoring equipment.
- Never open the battery cover during stimulation in order to avoid electrical shock.
- Turn off the stimulation or make sure that the amplitude for each channel is 0.0 mA before touching or removing the probes. Getting electrical stimulation through the fingers is unpleasant but not harmful.
- Observe caution when stimulating in the immediate vicinity of cellular phones that are switched on, since this may affect the output power of the stimulator.

3.1. CONTRAINDICATIONS

- **Implanted electronic devices.** Do not use the device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device, because this may cause electric shock, burns, electrical interference, or death.
- **TENS for Undiagnosed Pain.** Do not use the device as a TENS device on patients whose pain syndromes are undiagnosed.

3.2. WARNINGS

- Consult with physician. Consult with the patient's physician before using the device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Skin condition. Apply stimulation only to normal, intact, clean, healthy skin.
- Long term effects. The long-term effects of chronic electrical stimulation are unknown.
- Stimulation location. Stimulation over Neck or Mouth. Do not apply stimulation over the patient's neck (especially the carotid sinus) or the patient's mouth, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Stimulation across Chest. Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- Across the head. Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.
- Stimulation over compromised skin. Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- Stimulation near cancerous lesions. Do not apply stimulation over, or in proximity to, cancerous lesions.
- Stimulation over eyes. Do not apply stimulation directly on the eyes.
- Environment. Electronic monitoring equipment. Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Bath or Shower. Do not apply stimulation when the patient is in the bath or shower. Do not apply stimulation in humid atmosphere exceeding 75% of relative humidity. Sleeping. Do not apply stimulation while the patient is sleeping.
- Driving or operating machinery. Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
- Electrosurgical equipment or defibrillators. Disconnect the stimulation electrodes before using electrosurgical equipment or defibrillators. Otherwise skin burns may be caused below the electrodes and the device might be destroyed.
- Magnetic Resonance Imaging. Do not wear electrode or the device during Magnetic Resonance Imaging (MRI) scans as this may result in metal overheating and causing skin burns in the area of the electrode.
- Flammable or explosive environment. Do not use the device in areas where there is a risk of fire or explosion, such as oxygen-rich environments, in the vicinity of flammable anaesthetics, etc.
- Power supply. Never connect stimulation cables to an external power supply as there is a risk of electric shock.
- Near other equipment. Do not use the device beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the device works properly in the chosen configuration.
- Miscellaneous. Electrodes for Single Patient. Do not share electrodes with other persons. All users should have individual set of electrodes to prevent undesirable skin reactions or disease transmission.
- Accessories. Use this device only with the leads, electrodes, and accessories recommended by the manufacturer. Use of other accessories may adversely affect the performance of the device or may result in stronger electromagnetic emissions or reduce the electromagnetic immunity of the device.
- No Modification. No modification of the equipment is allowed.

3.3. PRECAUTIONS

- **Supervision.** Use this device only under the continued supervision of a licensed practitioner. Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.

- Abstract**

intend to use it for a prolonged period of time (more than 3 months). Extended storage of the batteries in the device might lead to batteries and device damage.

- Cable. The cable is best preserved if left attached to the stimulator between sessions. Do not jerk the cable or connection.
- Heat and Cold products. The use of heat or cold producing devices (e.g. electric heating blankets, heating pads or ice packs) may impair performance of the electrode or alter the patient's circulation/sensitivity and increase the risk of injury to the patient.
- Pulled muscles. Do not apply electrodes over pulled muscles. Using the stimulator on a previously extended muscle might further pull such muscle. The higher the stimulation intensity, the higher the risk to further overextends such muscle.
- Additional Precautions for TENS
 - TENS is not effective for pain of central origin, including headache.
 - TENS is not a substitute for pain medications and other pain management therapies.
 - TENS devices have no curative value.
 - TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Effectiveness of TENS is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.

3.4. DANGERS

- **Electrodes.** Any Electrode with a minimum active area of 16 cm² may be used with this device. Use of an electrode with an area less than 16 cm² can cause burns when the unit is used at higher intensities. Consult your clinician prior to using any electrode less than 16 cm².

3.5. ADVERSE REACTIONS

- Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin.
- Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.
- Patients may stop using the device and should consult with their physicians if they experience adverse reactions from the device.
- Precaution: Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first. Always turn off the stimulator before moving or removing any electrodes during a session.
- Precaution: Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/cm².
- Precaution: Do not apply stimulation in the vicinity of metal. Remove jewelry, body piercings, buckles or any other removable metallic product or device in the area of stimulation. Never use the electrodes contra-laterally, i.e. do not use two pins connected to the same channel on opposite segments of the body.
- Precaution: Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This reaction is connected with fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).
- Precaution: Sudden temperature changes can cause condensation to build up inside the stimulator. To prevent this, allow it to reach room temperature before use

2. INCREASE (left and right channel)

- ### 3. DECREASE (left and right channel)

- Decreases the amplitude (intensity of stimulation). Press and hold the button to decrease the amplitude continuously.
- Deactivates the key lock.
- Decreases the number of minutes when setting the timer (right button).
- Used when turning the program lock on and off (left button)*.
- Scrolls through choices in programming mode (left button)*.

4. PROGRAM

- Selects programs P1-P12.
- Pauses an ongoing program.
- Used when activating/deactivating the program lock*.

5. TIMER

- Initiates the timer setting.


6. PROGRAMMING/CONFIRMATION*

- Turns the stimulator into programming mode for custom programs P10–P12 when pressed for 2 seconds.
- Confirms settings in programming mode.

7. CABLE HOLDER

*Professional use

4. OVERVIEW/DISPLAY SYMBOLS

1. BROKEN CIRCUIT

Broken circuit. The reason for a broken circuit may be too high resistance or cable breakage. See chapter 10.

2. PAUSE

Paused program.

3. KEY LOCK

Activated key lock. The key lock is automatically activated if no key is pressed for 10 seconds. Deactivate the key lock by pressing the left or right **DECREASE** button.

4. PROGRAM LOCK

Activated program lock.

5. BATTERY STATUS

Empty batteries. This symbol is shown when the batteries are almost empty.

6. PROGRAM NUMBER

Selected program number.

7. REMAINING TIME

Remaining program time in minutes and seconds. Time is flashing during timer setting.

8. AMPLITUDE BARGRAPH (left and right channel)

Selected amplitude as a bargraph.

9. AMPLITUDE LEVEL (left and right channel)


Selected amplitude in mA.

10. STIMULATION/REST

Stimulation/rest indication for programs with intermittent stimulation. The upper part of the symbol is flashing during stimulation and the bottom part during rest time.

1. INSERT THE BATTERIES


A 1. Connect the electrodes to the cable(s).


A 2. P8-Pg: Attach the electrodes to your body


A 3. Connect the cable(s) to the CEFAR PERISTIM PRO.


B 1. Connect the probe(s) to the cable(s).

B 2. P1-P7: Insert the vaginal probe and/or the anal probe.

B 3. Connect the cable(s) to the CEFAR PERISTIM PRO.

TIMER

The preset treatment time is set to 30 minutes, but the timer function allows you to set the treatment time yourself. You have the option to change the timing from “- -” to 60 minutes. If you choose to stay on “- -”, the stimulation will go on until you stop the stimulation manually.

To set the timer:

1. Select a program and start the stimulation, see chapter OPERATION.
2. Press the **TIMER** button  to initiate the timer function.
3. Increase the time by pressing the right **INCREASE** button . The time is increased by 1 minute each time you press the button. Decrease the time by pressing the right **DECREASE** button . The time is decreased by 1 minute each time you press the button. The time is flashing on the display, while setting the timer.
4. Press the **TIMER** button  to confirm the timer setting.

PROGRAM PAUSE

You can pause the programs up to five minutes.

To pause a program:

1. If key lock is on, press any of the **DECREASE** buttons  to deactivate it.
2. Press the **PROGRAM** button  to pause the program.

Start stimulation again by pressing the **PROGRAM** button .


If stimulation is paused for more than five minutes, the stimulator turns off automatically to spare batteries.

6. REPLACEMENT OF BATTERIES

A battery symbol is shown on the display when the batteries are almost empty . As long as the stimulator is working normally you can continue the treatment. When stimulation feels weaker than usual or the stimulator turns off, it is time to replace the batteries.

If the stimulator is not used for some time (approximately 3 months), the batteries should be removed from the stimulator.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

REPLACING THE BATTERIES

1. Turn off the stimulator.
2. Locate the battery compartment on the back of the stimulator.
3. Remove the battery cover by sliding it downwards about 1 cm and then lifting it upwards.
4. Remove the batteries.
5. Insert the new batteries correctly according to the polarity markings (+ and -) inside the battery compartment, see picture.
6. Replace the battery cover.
7. Dispose of the exhausted batteries in accordance with local and national regulations.

NOTE! Non-rechargeable batteries may not be charged due to the risk of explosion.

7. PROGRAMS

PRESET PROGRAMS

CEFAR PERISTIM PRO offers eight preset programs for stress, urge, mixed and faecal incontinence treatment and one preset program for relief of bladder-related pain.

Programs 1-7 should be used with a vaginal and/or anal probe and programs 8-9 with surface electrodes.

In programs 1–6 intermittent stimulation is used, which means that stimulation time alternates with rest time without stimulation: stimulation/rest/stimulation. The rest time allows the muscles to recover.

Programs 1–3

STRESS INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P1	50 Hz	200 µs	3 s	6 s	30 minutes, 3–5 times a week
P2	50 Hz	200 µs	5 s	10 s	30 minutes, 3–5 times a week
P3	50 Hz	200 µs	10 s	20 s	30 minutes, 3–5 times a week

One cause of stress incontinence can be poorly functioning perineal muscles. Electrical stimulation of the perineal muscles with a vaginal and/or anal probe makes the muscles contract. The contractions help to identify and strengthen the weakened muscles. The stimulation must be as strong as possible without being painful, and it is helpful if the patient participates actively in the muscle contractions. The intensity should be sufficient enough to create a reflexive contraction of the anus. It is suggested that this training will be combined with the patient's own exercises for pelvic floor muscles.

Programs 4–6

MIXED INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P4	20 Hz	200 µs	3 s	6 s	30 minutes, 3-5 times a week
P5	20 Hz	200 µs	5 s	10 s	30 minutes, 3-5 times a week
P6	20 Hz	200 µs	10 s	20 s	30 minutes, 3-5 times a week

Program 4-6 is used with a vaginal and/or anal probe for treating mixed incontinence. The stimulation should be as strong as possible without being painful. It is helpful if the patient participates actively in the muscle contractions.

URGE INCONTINENCE

Programnummer	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P7	10 Hz	180 μ s	Continuous	–	30 minutes, 2–5 times a week

Normally the bladder muscle is relaxed between evacuations. In certain kinds of urge incontinence, however, the bladder is unstable and sudden uncontrolled bladder contractions occur during the "resting phase". Urge incontinence can be treated with low-frequency constant stimulation, which has a relaxing effect on the overactive bladder. Program 7 is intended for urge incontinence treatment and should be used with a vaginal and/or anal probe. The stimulation should be as strong as possible without being painful.

Program 1-6

FAECAL INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P1	50 Hz	200 µs	3 s	6 s	30 minutes, 3–5 times a week
P2	50 Hz	200 µs	5 s	10 s	30 minutes, 3–5 times a week
P3	50 Hz	200 µs	10 s	20 s	30 minutes, 3–5 times a week
P4	20 Hz	200 µs	3 s	6 s	30 minutes, 3–5 times a week
P5	20 Hz	200 µs	5 s	10 s	30 minutes, 3–5 times a week
P6	20 Hz	200 µs	10 s	20 s	30 minutes, 3–5 times a week

Faecal incontinence can be a result of e.g. weakened pelvic floor muscles, weakened or poorly functioning anal sphincter muscles or damage to the nerves controlling the anal sphincter muscles. In these cases electrical stimulation can be used with an anal probe to increase the tonus of the pelvic floor and the anal sphincter. The stimulation should be as strong as possible without being painful. The intensity should be sufficient enough to create a reflexive contraction of the anus. It is helpful if the patient participates actively in the muscle contractions.


Program 8

TENS (HIGH FREQUENCY) – PAIN RELIEF

Programnummer	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P8	80 Hz	180 µs	Continuous	-	At least 30 minutes daily

High-frequency TENS stimulation is an effective treatment for relief of both acute and chronic pain. External stimulation with surface electrodes can be used for diagnosed pain conditions in the lower part of the abdomen, as in cystitis.

The electrodes should be placed in the painful area over the pubic region, see electrode placement picture.


Program 9

TENS (BURST) – URINARY INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P9 Sacral nerves stimulation	2 Hz	180 µs	Continuous	–	At least 30 minutes daily
P9 Acupuncture pointsstimulation	2 Hz	180 µs	Continuous	–	30 minutes per session during at least one month. Start with 2–3 times a week and then reduce the treatment frequency.

Electrical stimulation for urinary incontinence is usually applied by a probe inserted in the vagina and/or anus. An alternative method is to use TENS to stimulate areas of the skin that are innervated by the same spinal cord segment (S2–S3) as the bladder and urethra. The stimulation is applied through surface electrodes placed on the skin between the anus and the genitals, alternatively in the lumbar region. In treatment of incontinence where the sacral nerves are stimulated, the electrodes should be placed in the region of S2–S3, see electrode placement picture. The stimulation should be strong enough to create a reflexive contraction of the anus.

Electrical stimulation of the posterior tibial nerve, via acupuncture points, has also proven to have a positive effect.

In treatment of incontinence where stimulation is applied over acupuncture points, the electrodes should be placed over SP6 and behind/beneath the medial malleolus, see electrode placement picture

CUSTOM PROGRAMS (P10-P12)

With the CEFAR PERISTIM PRO it is possible to create and store three custom programs for individual needs. To create a custom program, follow the programming procedure below. To use a custom program, follow the instructions in the chapter OPERATION. The program time is preset to 30 minutes but can be changed with the timer function, see section TIMER.

Programming

1. Press the **PROGRAM** button  repeatedly until Program 10, 11 or 12 is shown on the display. Select one of these programs.
2. Press the **PROGRAMMING/CONFIRMATION** button  for 2 seconds to enter the programming mode.
3. Press the **INCREASE (OR DECREASE)** button   to see the two different stimulation types alternate on the display: C, I. The programming chart below shows the setting possibilities for each stimulation type.

```

graph TD
    Root[Programming  
CEFAR PERISTIM PRO  
P10, P11, P12] --> C[Continuous  
(C)]
    Root --> I[Intermittent  
(I)]
    C --> F1[Frequency  
Preset: 80 Hz  
Range: 2-120 Hz]
    F1 -- "2-5 Hz" --> B[Burst]
    F1 -- "6-120 Hz" --> PW1[Pulse width  
Preset: 180 µs  
Range: 50-300 µs]
    I --> PW2[Pulse width  
Preset: 180 µs  
Range: 50-300 µs]
    PW2 --> F2[Frequency  
Preset: 10 Hz  
Range: 5-100 Hz]
    F2 --> ST[Stimulation time  
Preset: 5 s  
Range: 1-20 s]
    ST --> RT[Rest time  
Preset: 10 s  
Range: 1-30 s]
  
```

PROGRAM LOCK

The stimulator can be locked to prevent changing of programs.

To activate/deactivate the program lock:

1. Select the program you want to lock/unlock, see chapter OPERATION.
2. Press the **PROGRAM** button  and the left **DECREASE** button  simultaneously for 2 seconds.
3. Press the left **INCREASE** or **DECREASE** button  . "ON" is shown on the left side of the display when activating the program lock and "OFF" when deactivating it. (The button toggles between ON and OFF).
4. Press the **PROGRAM** button  to finish the program lock setting.


COMPLIANCE

Compliance gives you the possibility to monitor the use of the stimulator:

1. Turn the stimulator ON.
2. Press the **TIMER** button  and the right **DECREASE** button  simultaneously for 2 seconds.
3. The left side of the display shows the usage time in hours and the right side in minutes. To reset the usage time, press the right **DECREASE** button  for 2 seconds.
4. Wait for 5 seconds or press **TIMER** button .
5. The left side of the display shows the total usage time in hours and the right side in weeks. The total usage time cannot be reset.
6. Wait for 5 seconds or press the **TIMER** button  to exit compliance mode.

If you use the CEFAR PERISTIM PRO with surface electrodes the electrodes will eventually wear out and need to be replaced. It is recommended to replace the electrodes after approximately 20-40 times of usage.

For purchase information, contact your Cefar dealer or visit www.cefarcorp.com.

9. CARE INSTRUCTIONS

Taking care of and cleaning the Cefar-Compex equipment is simple with the following instructions:

- For care and cleaning instructions of probes, see probe package.
- Self-adhesive multi-use electrodes are re-moistened if necessary with a few drops of water and kept air-tight (in a plastic bag) on protective paper when they are not in use.
- Rinse the carbon rubber electrodes and the skin with water after use. Do not use detergent for the electrodes.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- Do not jerk cables or connections.
- The cables are best preserved if left attached to the stimulator between sessions.
- No other maintenance of the device is required. Service life of the device can vary depending on usage conditions. Typical service life is 7 years.
- Never service the device while in use
- The device should be operated in temperatures between 10°C and 40°C, atmospheric pressures between 50 and 106 kPa, and relative humidity between 30% and 75%.
- The device should be transported and stored in temperatures between -40 °C and 70 °C, atmospheric pressures between 50 and 106 kPa and relative humidity between 10% and 90%.
- Contact manufacturer for assistance in setting up, using or maintaining the equipment or report events.

- Check that all settings are correct (chapter OPERATION) and make sure that the electrodes are correctly placed.

- **Electrodes:** Slightly change the position of the electrodes.

- Check that the probe/electrode cable is properly connected to the stimulator.
- **Probes:** Make sure that the probe is correctly placed with good contact. In case of irritation/inflammation in the vagina or anus, please contact your health care provider.

- **Electrodes:** The skin is irritated. For advice on skin care, see chapter PRECAUTIONARY MEASURE.
- **Electrodes:** The electrodes begin to lose their stickiness and do not stick properly to the skin. Moisten the adhesive surface with a few drops of water before placing on the skin.
- **Electrodes:** The electrodes are worn out and need to be replaced.
- **Electrodes:** Slightly change the position of the electrodes.

- Check if the batteries need to be replaced, see chapter REPLACEMENT OF BATTERIES.
- **Electrodes:** Electrodes are too old and need replacement.

The broken circuit symbol indicates that the resistance is too high, or that a cable is broken.

- A too high resistance can be caused by a bad connection between the electrodes/probes and your skin/tissue, or that the electrodes need to be replaced.
- A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for the corresponding channel to 11 mA. If the amplitude now drops to 00.0 mA and $\mathcal{Q}_{\text{eq}}$ starts flashing, the cable needs to be replaced.

Note! Never increase the amplitude above 20 mA when you check for cable breaks, since this can damage the stimulator.

3

If the error symbol appears on the display when you start the stimulator, it means that the stimulator is broken and needs to be repaired.

Note! Do not use the stimulator – contact your Cefar dealer.

Cefar-Complex will only be responsible for service and repairs performed by Cefar or a distributor appointed by Cefar-Complex.

11. TECHNICAL DATA

CEFAR PERISTIM PRO is a stimulator with two non-independent channels intended for treatment of incontinence and bladder-related pain. It has nine preset programs and three custom programs.


Treatment with electrical stimulation requires the stimulation current to penetrate the resistance of the skin and the electrode, about 1000 ohms. CEFAR PERISTIM PRO can penetrate this resistance and maintain a current of up to 99.5 mA. With a change in load from 100 to 1000 ohms, the stimulation current changes less than 10% from the set value.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

CEFAR PERISTIM PRO

Number of channels	2 (non-independent)
Constant current	Up to a resistance of 1000 ohm (increased load can reduce the maximum current)
Stimulation current/channel	0–99.5 mA (charge maximum: 30µC)
Waveform	Symmetrical biphasic pulse, 100% compensated
Number of preset programs	9
Number of custom programs	3
Stimulation forms	Continuous stimulation Intermittent stimulation
Max pulse duration	300 µs
Max frequency	120 Hz
Timer	1 to 60 min/Off
Environment for storage, use and shipping	Temperature 10° C–40° C Air humidity 30%–75% Air pressure 700 hPa–1060 hPa
Power source	2 x 1.5 V AA non-rechargeable or 2 x 1.2 V AA rechargeable batteries
Current consumption for one channel, 80 Hz, 30 mA	150 mA
I r.m.s. max/channel	27 mA
Size	120 x 50 x 30 mm
Weight	Ca. 180 g

12. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should ensure that it is used in such an environment		
Emissions Tests	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The CEFAR PERISTIM PRO Electrotherapy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	CEFAR PERISTIM PRO Electrotherapy System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not Applicable - Battery powered	
Voltage fluctuations IEC 61000-3-3	Not Applicable - Battery powered	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment specified below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Risk assessment on the CEFAR PERISTIM PRO Electrotherapy System indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not Applicable - Battery powered Not Applicable - signal lines less than 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 1kV differential mode (line to line) + 2kV common mode (line to ground)	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR PERISTIM PRO Electrotherapy System requires continued operation during power mains interruptions, it is recommended that the CEFAR PERISTIM PRO Electrotherapy System be powered from an uninterrupted power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the AC mains voltage prior to application of the test level.

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	[V ₁] V, where V ₁ = 3V	Portable and mobile RF communications equipment should be used no closer to any part of the CEFAR PERISTIM PRO Electrotherapy System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{35}{V_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{35}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CEFAR PERISTIM PRO Electrotherapy System is used exceeds the applicable RF compliance level above, the CEFAR PERISTIM PRO Electrotherapy System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CEFAR PERISTIM PRO Electrotherapy System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V] V/m.

Rated maximum output power of transmitter P (W)	Separation distance according to frequency of transmitter d (m)		
	150 kHz to 80 MHz $d = \frac{35}{V_1} \sqrt{P}$ (where $V_1 = 3V$)	80 MHz to 800 MHz $d = \frac{35}{E_1} \sqrt{P}$ (where $E_1 = 3V/m$)	800 MHz to 2,5 GHz $d = \frac{7}{E_1} \sqrt{P}$ (where $E_1 = 3V/m$)
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Part 15 of the FCC Requirements

- This device may not cause harmful interference
- This device must accept any interference received, including the interference that may cause undesired operation.

T9J-RN42

6514A-RN42

13. REFERENCES

Fall M, Lindström S

Electrical stimulation: A physiologic approach to the treatment of urinary incontinence.

Journal of Urologic Clinics of North America; 2, Vol 18: 393-407, 1991

Fall et al. Transcutaneous Electrical Nerve Stimulation in classic an nonulcer interstitian cystitis.

Urologic clinics of North America Vol 21, No 1, Feb 1994

Walsh et al. Non-invasive antidromic neurostimulation, a simple effective method for improving bladder storage. *Neurourol Urodyn* 2001;20(1):73-84

Amarengo et al. Urodynamic effect of acute transcutaneous posterior tibial nerve stimulation in over-active bladder. *J Urol.*2003 Jun;169(6):2210-5

Zöllner-Nielsen M, Samuelsson S.M

Maximal electrical stimulation of patients with frequency, urgency and urge incontinence. Report of 38 cases. *Acta Obstet Gynecol Scand* 1992;71:629-631

Sand P.K, Richardson D.A Et al.

Pelvic floor electrical stimulation in the treatment of genuine stress incontinence: A multicenter, placebo-controlled trial. *Am J Obstet Gynecol* Vol 173, No 1, 1995

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© 2016 DJO - 1566 - Rev A