Gebrauchsanleitung (DE-1 - DE-27)  
Instructions for Use (EN-28 - EN-53)  
Mode d'emploi (FR-54 - FR-79)  
Istruzioni per l'uso (IT-80 - IT-107)

# Domotherm®

**UEBE**  
Germany  
[www.uebe.com](http://www.uebe.com)

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## Table of contents

<b>A Safety instructions</b>	30
<b>B Operating the unit</b>	
1. Advantages of the infrared ear thermometer	32
2. Fever classification	33
3. Description of the unit	34
4. Control displays and symbols	35
5. Start-up of the unit	36
6. Battery change	36
7. Temperature measurement	37
8. To avoid inaccurate measurements	39
9. Call up measured values	40
10. Conversion from Fahrenheit to Celsius degrees and vice versa	41
<b>C Technical information</b>	
1. Error messages	42
2. Customer service	43
3. Technical data	44
4. Metrological control	45
5. Applied standards	45
<b>D General provisions</b>	
1. Explanation of symbols	46
2. Disposal	46
3. Electromagnetic compatibility (EMC)	47
<b>E Maintaining the unit</b>	51
<b>F Warranty</b>	52

## Introduction

The Digital Infrared Ear Thermometer Domotherm® Ear (hereinafter also referred to as unit) is a high-quality product providing the latest achievements of technology and tested in accordance with the international standards.

By its unique technology, the Domotherm® Ear facilitates solid measured values free of any disturbing thermal influence with every measurement. Each time it is switched on, the unit performs an auto-test to guarantee the indicated measuring accuracy.

As with any other thermometer, the right measuring technique is an essential requirement for exact measurement results. Please read the operating manual carefully before using the unit and keep it in a safe place.

### Why measure the temperature in the ear?

Temperature measurement is about measuring the body temperature that corresponds to the temperature of the internal organs.

A temperature measured in the ear reflects the body's core temperature very precisely because the eardrum and the temperature control centre in the brain, the hypothalamus, are supplied by joint blood vessels.

Thus, changes in the body temperature are indicated faster and more precisely in the ear than in other places.

### Mode of operation:

The unit measures the infrared energy emitting from the middle ear and the surrounding tissue. This energy is absorbed by lenses and converted into temperature values.

The measured values directly received from the middle ear (eardrum) guarantee the most precise ear temperature.

## A Safety instructions

### Important safety instructions

- Please use the thermometer only for the purposes it is intended for. The unit is only suitable for temperature determination in the acoustic meatus. Please do not use it for rectal, oral and axillary temperature measurement. Please note the general safety precautions when applied on children.
- Do not immerse this thermometer in water or other liquids (not waterproof). With respect to cleaning and disinfection, please follow the instructions in the chapter «Maintaining the unit».
- Do not expose the unit to direct sunlight and keep it in a dust-free and dry place. Please note the operation, storage and transport temperature.
- Do not use the thermometer if indications of damages are visible on the measuring part or the unit itself. Should the thermometer really be damaged: Please do not try to repair it yourself! Please contact your UEBE Medical GmbH customer service in this case.
- Ear wax in the acoustic meatus and on the measuring probe can influence the temperature determination. Please make sure that no impurities can affect the measurement result prior to the measurement.
- This thermometer consists of high-quality precision parts. Do not drop the unit! Protect it from severe impacts and shocks or other unfavourable impairments. Please also refer to storage and transport conditions.
- Do not use the thermometer in an inflamed or injured ear.
- In case of storage outside the ambient operating temperature, the thermometer must be stored unpacked in the treatment room for at least 4 hours to equal the temperature. We recommend to customers storage in living areas for immediate use.

EN-30

## Safety instructions

### A

#### Caution:

- Use of this infrared ear thermometer does not replace appointment with a doctor.
- The thermometer is not waterproof! Please do NOT immerse in liquids!
- The thermometer is not suitable to determine the basal body temperature.

EN-31

## B Operating the unit

### 1. Advantages of the infrared ear thermometer

- Free of protective cover: It is not necessary to use a protective cover for temperature measurement which makes the thermometer user-friendlier and more cost-effective.
- Quick measurement: the innovative infrared technology allows for ear temperature measurement within only 1 second.
- Precise and reliable: The unique construction of the measuring probe and of the infrared sensor and the calibration methods facilitate very precise and reliable ear temperature measurements with this unit.
- Simple and convenient use: Due to its ergonomic special design, the thermometer can be used in a simple and convenient manner. Measurements can be performed even if the child is asleep. This thermometer turns out to be especially convenient for use with children as it is more comfortable than the rectal or oral measurement.
- Measured value memory: The unit is very user-friendly and automatically shows the value of the last measurement when it is activated.
- Display of multiple measured values: The user can call up the last 12 measured values stored so that they are able to better keep track of the temperature fluctuations.

## Operating the unit

B

- Safe and hygienic
  - No risk for breakage of glass or mercury ingestion.
  - Absolutely safe use with children.
  - Clean the measuring probe with one of the cleaning pads attached or with a paper tissue slightly soaked in alcohol which makes the thermometer 100% hygienic.
- Fever alarm: 10 short acoustic signals indicate that the temperature of the patient is more than 37.5°C.
- The thermometer may not be used with a protective cover to determine the body temperature because this would falsify the measurement results!

### 2. Fever classification

- 37.5°C to 38.0°C elevated temperature
- 38.0°C to 39.0°C moderate fever
- as of 39.0°C high fever

## B Operating the unit

### 3. Description of the unit

- ① Protective cap
- ② Measuring probe
- ③ LCD display
- ④ O/I button (ON/OFF)
- ⑤ Battery cover
- ⑥ Start button

EN-34

## Operating the unit

B

### 4. Control displays and symbols

LCD display	Meaning of the display	Description
	Display test: Indication of all display elements	Switch the unit on with the O/I button and all display segments are shown for 2 seconds.
	Memory	The value of the last measurement appears automatically for about 2 seconds.
	Ready	The unit is ready for measurement and the °C or the °F symbol, as the case may be, flashes.
	Measuring process completed	The measured value appears on the LCD display with illuminating °C or °F symbol, as the case may be. The flashing temperature symbol indicates readiness for the next measurement. Additionally, the background is illuminated for about 5 seconds.
	Indication of a temperature outside the ear	The symbol of a crossed out ear appears on the LCD display if the measured value is outside the temperature range of 32.0 to 42.2°C (89.6 to 108°F).
	Battery change	When the unit is activated, the battery symbol appears continuously. Measurement is no longer possible, the battery must be changed.

## B Operating the unit

### 5. Start-up of the unit

- Open battery case
- Remove contact protecting strip
- Close battery case

### 6. Battery change

- Open battery case
- Push old battery with a thin item or screwdriver in the empty space provided forward in the direction of the start button.. Remove battery in upward direction. Place the new battery (+pole facing upwards) under the metal brackets. Press the battery slightly in direction of the start button direction until it clicks into place under the two plastic brackets.
- Close battery case.
- Please note:

#### Change the batteries

- if the empty battery symbol appears in the display after the display test.
- if nothing is indicated in the display despite of pressing the O/I button (ON/OFF).

#### Note:


The battery sign is also shown during the display test after activation, then the sign disappears if the battery is alright.

## Operating the unit


### B

### 7. Temperature measurement

- Press the O/I button. The display test is activated.

- The value of the last measurement appears automatically with the symbol «M» (Memory) on the display.
- Before the °C or the °F symbol flashes, there is a beep by which the thermometer indicates its readiness for measurement.
- Straighten the acoustic meatus by pulling the ear upwards and backwards so that the middle ear is clearly visible.
  - Children under 1 year:  
Pull the ear straight backwards.

- Children from 1 year and adults:  
Pull the ear backwards and upwards.

## B Operating the unit

- While you are holding the ear in its position, insert the measuring probe into the acoustic meatus. Press the start button when the thermometer is placed on the point of measurement.

Keep pressing the button as long as the «beep» signal appears. This signal confirms completion of the measurement.

### Important note:

Press the start button only if the thermometer is finally placed on the point of measurement. Otherwise the taking of the core temperature could be missed due to the quick measuring process.

- Pull the thermometer out of the acoustic meatus. The LCD display indicates the measured temperature.

### Note:

- As of 37.5°C (99.5°F) it is indicated with 10 short acoustic signals that there is elevated temperature. (Fever alarm)
- In order to achieve highest accuracy in subsequent measurements, please wait for 30 seconds after each measurement.
- If ear wax accumulates on the measuring tip this can lead to inaccuracies of the temperature measurement or to a cross infection between different users. For this reason, please provide for a clean measuring probe after each measurement. Please follow the instructions in the chapter «Maintaining the unit».
- After cleaning the measuring tip, please wait for 5 minutes before the next use so that the thermometer achieves its operating temperature.

## Operating the unit

- It is best to lay toddlers flat on their stomach with their head to one side so that the ear shows in upward direction. With older children or adults it is better to stand slightly sideward behind the patient.
- Always measure the temperature in the same ear as the temperature values can vary from ear to ear.
- After sleep, please wait a couple of minutes before you measure the ear temperature.
- In the following situations it is recommended to perform three temperature measurements in the same ear and to take the highest measured values of them:
  - Newborns in the first 100 days.
  - Children under three years or patients with endangered immune system, in which case the existence or non-existence of fever can be decisive.
  - Users should practice the correct use of the infrared ear thermometer until they get acquainted with the unit and receive consistent measured values.

### 8. To avoid incorrect measurements:

Measurements taken on surrounding tissue of the acoustic meatus result in lower measured values and can cause an incorrect fever diagnosis.

Please refer to the chapter „Temperature measurement”.

This thermometer was clinically tested and turned out to be safe and precise if used according to the operating manual.

## B Operating the unit

### 9. Call up measured values

This thermometer allows to store and call up the last 12 measured values. To do so, please proceed as follows:

LCD display	Meaning of the indication	Description
	Memory call	Press the start button for about 1 second, release button, the symbol «M» flashes.
	Call 1 (the measured value last identified)	Press the start button again and release to call up the value last measured. On the display appears «1» jointly with the memory symbol.
	Call 2 (the second last measured value)	Press the start button repeatedly and release to call up the second last measured value.
	Call 12 (the oldest of 12 measured values)	Press the start button and release several times consecutively to call up the last 12 values saved.

When you press the start button again after 12 called up measured values and release, the course of this sequence starts again, i.e. with measured value 1.

## Operating the unit

### B

### 10. Conversion from Fahrenheit to Celsius degrees and vice versa

This thermometer can indicate the measured temperature either in Fahrenheit or in Celsius degrees.

To convert the display from °C to °F, simply switch off the unit (O/I) and press the start button for about 5 seconds until three lines and °C appear in the display.

By pressing the start button again, you can now switch between °C and °F.

After switching off the unit, the setting last selected remains active.

## C Technical information

### 1. Error messages

LCD display/problem	Meaning of the indication	Possible cause and error removal
	Measurement temperature too high	Display «H» if the measurement temperature is more than 100.0°C or 199.9°F.
	Measurement temperature too low	Display «L» if the measurement temperature is less than 0°C or 32.0°F.
	Ambient temperature too high	Display «H» in connection with «▲» if the ambient temperature is more than 40.0°C or 104.0°F. Consider the item „Ambient operating conditions“ in the field Technical data.
	Ambient temperature too low	Display «L» in connection with «▼» if the ambient temperature is less than 5°C or 41.0°F. Consider the item „Ambient operating conditions“ in the field Technical data.
	Malfunction indicator	Switch off unit, remove battery for 2 minutes, then make a new attempt. If same problem appears, call customer service.

## Technical information

### C

	„Empty“ indicator	Please check if the battery is positioned correctly. Also check the polarity (<+> and <->) of the battery. See chapter Battery change.
	Battery weak	If the battery sign flashes during the measurement process, the battery is weak; please exchange batteries.
	Battery empty indicator	If the battery symbol appears continuously and as a single symbol, the battery must be exchanged as soon as possible. See chapter Battery display and battery change.

### 2. Customer service

Only the manufacturer or an expressly authorised body may repair the unit.

Please contact:

UEBE Medical GmbH  
Zum Ottersberg 9  
97877 Wertheim, Germany  
Phone no.: +49 (0) 9342/924040  
Fax no.: +49 (0) 9342/924080  
E-mail: info@uebe.com  
Internet: www.uebe.com

## C Technical information

### 3. Technical data

Type:	Infrared ear thermometer free of protective cover
Measuring range:	0°C to 100.0°C (32.0°F to 212.0°F)
Measuring accuracy:	Medical sector: 35.0°C to 42.0°C ±0.2°C 95°F to 107.6°F ±0.4°F
Measuring accuracy:	Outside the medical sector: 0°C to 32°C as well as 42°C to 100°C ±1.0°C 32°F to 89.6°F as well as 107.6°F to 212°F ±33.8°F 32°C to 34.9°C as well as 42°C to 42.2°C ±0.3°C 89.6°F to 94.8°F as well as 107.6°F to 108°F ±0.6°F
Display:	Liquid crystal display with display unit 0.1°C (0.1°F)
Beep signals:	a. 1 short acoustic signal: The unit is switched ON and ready for measurement. b. 1 long acoustic signal (1 s.): Measuring process completed if the measured value is under 37.5°C (99.4°F). 10 short acoustic signals: If the measured value is higher than or equals 37.5°C (99.4°F). c. 3 short acoustic signals: System error or malfunction.
Memory:	a. Value last measured appears as soon as unit is switched on. b. 12 measurements by memory call
Ambient operating conditions:	6°C to 40°C (41.0°F to 104°F) Rel. air humidity 15% - 95%
Storage and transport conditions:	-25°C to 55°C (-13°F to 131°F) Rel. air humidity 15% - 95%
Automatic switch-off:	About 1 minute after the last measurement. The unit switches off automatically to protect against misuse.

## Technical information

### C

Serial number:	On the unit there is a serial number <u>SN</u> that clearly identifies the unit.
Battery:	CR2033
Size:	153mm (L) x 31mm (W) x 40mm (H)
Weight:	53g (with battery), 50g (without battery)

### 4. Metrological control

A metrological control is generally recommended every 2 years. However, in Germany, expert users are required to do so according to the „Medizinprodukte-Betreiberverordnung“ [Medical devices operator ordinance].

It can be performed either by the UEBE Medical GmbH, by an authority responsible for metrology or by authorised maintenance services. Please refer to your national requirements.

**Important:** Without the manufacturer's consent, no changes may be made on this unit (e.g. opening of the unit other than for battery change purposes).

### 5. Applied standards

- Meets the requirements ASTM E-1966
- EN 12470:2003 Medical thermometers part 5: Requirements to infrared ear thermometers
- IEC 60601-1:2005 ·Medical electrical devices part 1: General determination for safety including essential performance features
- IEC 60601-1-2 : 2001 + A1:2004 ·Medical electrical devices - part 1-2: General determination for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests
- This unit corresponds to the Council Directive 93/42/EEC from 14 June 1993 about medical products and bears the sign CE 0123 (TÜV SÜD Product Service GmbH)

## D General provisions

### 1. Explanation of symbols

**CE0123** This unit corresponds to the Council Directive 93/42/EEC from 14 June 1993 about medical products and bears the sign CE 0123 (TÜV SÜD Product Service GmbH)  
Devices with CE marking are quality-controlled according to this directive and provide higher accuracy than the former calibration.

Degree of protection against electric shock: TYPE BF

Please note operating manual

55°C Storage and transport conditions ambient temperature  
-25 to 55°C (-13°F to 131°F)

95% Protect against moisture/humidity  
Rel. air humidity 15% - 95%

Keep dry

Manufacturer

### 2. Disposal

Batteries and technical devices are not household waste and must be delivered to the corresponding collection and disposal sites.

EN-46

## General provisions

### D

### 3. Electromagnetic compatibility (EMC) Technical description

The unit complies with the EMC requirements of the international standard EC60601-1-2. The requirements are fulfilled under the conditions described in the following.

The unit is an electric medical product and is subject to special precautions with respect to EMC that must be published in the operating manual.

Portable and mobile HF communication facilities can affect the unit. Use of the unit with accessories not approved can affect the unit in a negative way and change the electromagnetic compatibility.

The unit should not be used immediately next to or between other electric devices.

Table 201 - Guidance and manufacturer's declaration - electromagnetic emissions

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	N/A	

EN-47

## D General provisions

Tabelle 202 – Guidance and manufacturer's declaration - electromagnetic immunity

Electromagnetic environment - guidance	Compliance level	IEC 60601 test level	Immunity test
Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Electrostatic discharge (ESD) IEC 61000-4-2
N/A	N/A	±2 kV for power supply lines ±1 kV for input/output lines	Electrical fast transient/burst IEC 61000-4-4
N/A	N/A	±1 kV differential mode ±2 kV common mode	Surge IEC 61000-4-5
N/A	N/A	<5% U <sub>r</sub> (>95% dip in U <sub>r</sub> ) for 0,5 cycle 40% U <sub>r</sub> (60% dip in U <sub>r</sub> ) for 5 cycles 70% U <sub>r</sub> (30% dip in U <sub>r</sub> ) for 25 cycles <5% U <sub>r</sub> (>95% dip in U <sub>r</sub> ) for 5 sec	Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11
Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.	3 A/m	3 A/m	Power frequency (50/60 Hz) magnetic field IEC 61000-4-8
NOTE U <sub>r</sub> is the a.c. mains voltage prior to application of the test level.			

## General provisions

## D

Table 204 - Guidance and manufacturer's declaration - electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of Domotherm® Ear, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance N/A d=1,2/P 80 MHz to 800 MHz d=2,3/P 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	b Interference may occur in the vicinity of equipment marked with the following symbol: 
			NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radios broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Domotherm® Ear is used exceeds the applicable RF compliance level above, Domotherm® Ear should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Domotherm® Ear			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.			

## D General provisions

Table 206 - Recommended separation distances between portable and mobile RF communications equipment and Domotherm® Ear

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz N/A	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2,5 GHz $d=2.3\sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using

the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.


NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## E Maintaining the unit

E

Please use one of the alcohol tissues provided or a piece of cotton wool moistened with alcohol (70% isopropyl) to clean the thermometer casing and the measuring probe. Make sure that no liquid gets into the inside of the thermometer. Do not use scratching detergents, dilutors or benzol to clean and do not immerse the unit in water or other cleaning liquids.

Take care not to scratch that the LCD surface. Remove the battery from the unit if you do not use it for a longer period of time so that the thermometer is not damaged by a leaking battery.

## F Warranty

### Warranty

The unit was manufactured and tested with care. In case it nevertheless shows defects upon delivery, we give a warranty according to the following conditions:

- During the warranty period of 2 years as of date of purchase, we remove such defects at our discretion and at our costs by repair (after return to our factory) or replacement delivery of a defect-free unit.
- Normal wear and tear of wear and tear elements as well as damages caused by non-adherence to the operating manual, improper use (e.g. unsuitable power sources, breakages, leaking batteries) and/or disassembly of the unit by the customer are not covered by the warranty.  
Furthermore, the warranty does not establish claims for damages against us.
- Warranty claims can only be asserted during the warranty period and by presentation of the sales receipt. In case of a warranty claim, please send the unit jointly with the sales receipt and the description of the complaint to:

UEBE Medical GmbH  
Zum Ottersberg 9  
97877 Wertheim  
GERMANY

- The statutory claims of the purchaser against the seller based on defects in accordance with §437 BGB [German Civil Code] are not limited.

Detailed user information about our Domotherm® thermometer, visomat blood pressure measurement units and Cyclotest products as well as services is provided on our website [www.uebe.com](http://www.uebe.com).

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(C) Copyright 2009 UEBE Medical GmbH