

INTELECT® RPW LITE

EASY ACCESS TO RADIAL PRESSURE WAVE THERAPY

EFFICIENCY

Treatment outcomes are determined by the level of energy released by the high quality transmitters. The Intelect® RPW Lite is delivered with the D20-S transmitter. D20-T, Ro40, C15 and D115 are available as an option.

PORTABILITY

Compact and easy to handle.

INDICATIONS

- ACHILLES TENDINOPATHY (1)
- CARPAL TUNNEL SYNDROME (2)
- PLANTAR FASCIITIS (3)
- TRIGGER FINGER (4)
- TRIGGER POINTS (5)
- PAIN AND INFLAMMATION MANAGEMENT IN ORTHOPAEDIC CONDITIONS (6)

TECHNICAL INFORMATION

30680-230V: 230VAC 50-60Hz 700VA

30680-115V: 115VAC 60Hz 700VA

30680-100V: 100VAC 60Hz 700VA

DIMENSIONS: 289x238x310mm

WEIGHT: 9.8Kg

ENERGY: 6 Energy levels: 6-18Hz/max 2.7 bar

CONFORMITY: this device complies with the applicable standard: FN/IFC 60601-1 CAN/CSA-C22 2 No 6011 UL Std. No 60601-1

CERTIFICATION: CE 0197

CLASSIFICATION ACCORDING TO MDD: Class IIa device 93/42 60601-1



- 1. Gerdesmeyer L. Current evidence of extracorporeal shock wave therapy in chronic Achilles tendinopathy. Int J Surg. 2015 "
- 2. Wu YT. Effect of radial shock wave therapy for carpal tunnel syndrome: .J Orthop Res. 2016 Jun
- $\textbf{3.}\,\textbf{Ogden JA.}\,\textbf{Shockwave the rapy for chronic proximal plantar fasciitis Foot Ankle Int.}\,\textbf{2002}\,\textbf{Apr}$
- $4. Yildirim\ P. Extracorporeal\ shock\ wave\ the rapy\ versus\ corticosteroid\ injection\ in\ the\ treatment\ of\ trigger\ finger:\ . J\ Hand\ Surg\ Eur\ Vol.\ 2016\ Nov\ P. F.\ And\ Surg\ Eur\ Vol.\ 2016\ Nov\ P.\ And\ Surg\ Eur\ P.\ And\ Sur$
- 5. Nedelka T.Mechano-transduction effect of shockwaves in the treatment of lumbar facet joint pain:.Neuro Endocrinol Lett. 2014
- 6. Schmitz C. Efficacy and safety of extracorporeal shock wave therapy for orthopedic conditions.Br Med Bull. 2015



INTELECT® RPW LITE

ORDERING INFORMATION

PART NUMBER	DESCRIPTION
30680-230V	Intelect® RPW Lite 230V set including standard accessories
30680-115V	Intelect RPW Lite 115V unit set including standard accessories
30680-100V	Intelect RPW Lite 100V unit set including standard accessories

STANDARD ACCESSORIES

PART NUMBER	DESCRIPTION
30885-1	Intelect RPW Lite SPARROW™ handpiece kit grey, including 1x projectile; 1x guiding tube; 1x wrench; 1x allen key
30873	Transmitter D20-S Black RPW Lite, 20mm
-	1x Power cord
30891	Operating manual and quick start guide
22601	Transmission gel, 250 ml
28739	Guiding tube cleaning brush
30890	Tube set (includes 3 tubes)

OPTIONAL ACCESSORIES

PART NUMBER	DESCRIPTION
30885-1	Intelect RPW Lite SPARROW handpiece kit grey, including 1x projectile; 1x guiding tube; 1x wrench; 1x allen key
30885-2	Intelect RPW Lite SPARROW handpiece kit green, including 1x projectile; 1x guiding tube; 1x wrench; 1x allen key
29104	SPARROW handpiece revision kit, including 2x projectiles; 2x guiding tubes short; 2x o-ring set; 1x o-ring guide
30873	Transmitter D20-S Black RPW Lite, 20mm
30876	Transmitter D20-T Black RPW Lite, 20mm
30842	Transmitter Ro40 Black RPW Lite, 15mm
30875	Transmitter C15 Black RPW Lite, 15mm
30874	Transmitter DI15 Black RPW Lite, Deep Impact 15mm
30883	Handpiece holder black
30890	Tube set (includes 3 tubes)
22601	Transmission gel, 250 ml
28739	Guiding tube cleaning brush
0.0032.012	Power cord EU
0.0032.012-US	Power cord USA
13-7458-CN	Power cord China/Australia/New Zealand
13448-CH	Power cord Switzerland



DJO® Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM