

INTELECT® **RPW 2** GENERATION 2020

0

Chattanooga®

INTELECT[®] **RPW 2** GENERATION 2020

With 10 years of experience and 10,000 RPW devices sold all over the world Chattanooga® has learned a lot about what really matters to you and your patients and integrated it into our next generation clinical Radial Pressure Wave unit.

The new Intelect® RPW 2 combines streamlined design with ergonomics and patient-comfort with ease-of-use resulting into the most compact clinical unit currently on the market.

With "made in Germany & the USA" as a quality label we are proud to offer you the latest technology in a device you can trust!

KEY FEATURES



SUGGESTED **PROTOCOLS LIBRARY**

Available at a finger click; treatment pictures, anatomical and pathological library, treatments guidelines.

LARGE GRAPHIC COLOR TOUCH **SCREEN**

High resolution 7" capacitive touch screen for greater customer experience.

tanoog

NEW FALCON® HANDPIECE

allowing lower Bar pressures for an increased range of indications and applicators.

COMFORT MODE

Unique automatic energy ramp up designed to provide better patient comfort.

V-ACTOR® HF*

Vibration frequency up to 50Hz for Vibration therapy and muscle relaxation.





LARGE CHOICE OF TRANSMITTERS*

17 different transmitters available for tendinopathies, trigger points, spine and fascia treatments.

EASY FIT FOR EVERY TREATMENT ROOM

Slim design 41x42x110cm (LxWxH), not a lot of space needed for this new RPW generation.



ERGONOMIC DESIGN

What you is need close in hand; handpiece, transmitter tray, gel bottle, cables holders, soft pouch.

EXTENDED PRESSURE RANGE

The Falcon handpiece provides a greater pressure range from 0.3 to 5 bar, making it more comfortable for the patient and allowing Fascia and Spine treatments. NOTE: PERI-V-ACTOR® and SPINE-V-ACTOR® can only be used under 3 Bar pressure.

* Optional accessories

TARGETED PAIN RELIEF

The Intelect® RPW 2 allows for the treatment of indications with radial, pneumatically-generated, low-energy acoustic waves, or 'pressure pulses'. The applied pressure pulse propagates radially within the tissue, generating a therapeutic effect.

The V-ACTOR® HF optional hand piece allows mechanical vibration, for vibration massage

INDICATIONS

- ACHILLES TENDINOPATHY (1)
- DISORDERS OF TENDONS INSERTIONS (2)
- MYOFASCIAL TRIGGER POINTS (3)
- PAIN & INFLAMMATION IN ORTHOPAEDIC CONDITIONS (4)
- PLANTAR FASCIITIS (5)
- PULSE VIBRATION MASSAGE (6)

RECOMMENDED TREATMENT GUIDELINES

- LOCALIZE THE PAINFUL AREAS YOU PLAN TO TREAT
- SELECT MOST APPROPRIATE TRANSMITTERS BASED ON TREATMENT GOAL AND TARGETED TISSUE
- APPLY ULTRA SOUND WATER BASED GEL ON THE SKIN OVERLYING THE TREATMENT AREA
- ALWAYS START WITH LOWER LEVEL OF INTENSITY
- PERFORM A POST-TREATMENT TEST TO SEE FOR BETTER FUNCTION AND DECREASED IN PAIN

As a result of clinical evidence, most treatment protocols are based on 5 to 10 sessions of 1,500-2,500 shocks

EDUCATION

Access to an online education platform and scientific research library is offered to every Chattanooga® Intelect® RPW 2 user

CLINICAL INFORMATION

 Gerdesmeyer L. Current evidence of extracorporeal shock wave therapy in chronic Achilles tendinopathy.NCBINCBI Logo Int J Surg. 2015 Dec;24(Pt B):154-9. doi: 10.1016/j.ijsu.2015.07.718. Epub 2015 Aug 29."

Gerdesmeyer L. Current evidence of extracorporeal shock wave therapy in chronic Achilles tendinopathy.NCBINCBI Logo Int J Surg. 2015

3. Nedelka T.Mechano-transduction effect of shockwaves in the treatment of lumbar facet joint pain.Neuro Endocrinol Lett. 2014

 Schmitz C. Efficacy and safety of extracorporeal shock wave therapy for orthopedic conditions: a systematic review on studies listed in the PEDro database. Br Med Bull. 2015;116:115-38. doi: 10.1093/bmb/ldv047. Epub 2015 Nov 18.

max. 3 bar

5. Ogden JA. Shockwave therapy for chronic proximal plantar fasciitis: a meta-analysis.Foot Ankle Int. 2002 Apr;23(4):301-8

6. S. Broadbent, JJ. Rousseau, RM. Thorp, SL. Choate, FS. Jackson, DS. Rowlands. Vibration therapy reduces plasma IL6 and muscle soreness after downhill running 2010 Sep;44(12):888-94

INTELECT® RPW 2 TRANSMITTERS

A choice of standard transmitters for Radial Pressure Wave (RPW) indications. Each transmitter is unique by shape and material to deliver it's own energy flux density and wave profile.

TRIGGER POINTS / TENDINOPATHIES



15 mm Energy beam transmitter with concave coupling surface, best for pain zones near the skin surface Penetration depth 0-35mm Intensity level: Medium



R15 15mm transmitter for any type of tendinopathies Penetration depth 0-35mm Intensity level: Medium



C15 CERAma-x[®] Ceramic Energy 15mm transmitter for any type of tendonopathies Penetration depth 0-35mm Intensity level: High



Golden Depth 15mm Deep Impact®transmitter for deep target areas, chronic disorders, local trigger points Penetration depth 0-50mm Intensity level: High



A6 Acupuncture like transmitter 6mm for local trigger points Penetration depth 0-20mm Intensity level: Low



T10 Finger transmitter 10mm for deep target tendonopathies Penetration depth 0-40mm Intensity level: High



F15 White soft transmitter 15mm for superficial pain regions, muscles of mastication, cervical spine Penetration depth 0-30mm



ATLAS Soft flex

Soft flexible tipped transmitter 15mm for highly sensitive regions and triggerpoints in the cervical spine Penetration depth 0-15mm Intensity level: Low

Standard Oscillator, 20mm transmitter or muscle and connective

MYOFASCIAL TREATMENTS

Intensity level: Very Low







D2O-T Golden Oscillator 20mm transmitter for professional athletes, High energy application Penetration depth 0-50mm Intensity level: Very High



D35-S

Large oscillator 35mm for large muscles and connective tissue Penetration depth 0-60mm Intensity level: Medium

SPINE TREATMENTS

Double Tip shape to allows treatments of the paravertebral region

SPINE-ACTORS[®] 3 sizes span from 20 to 35mm



FASCIA TREATMENTS

Unique shape especially designed to dissolve clots in the Fascia structures



PERI-ACTOR® KNUCKLE

For Medium-sized and large fasciae thoracolombar fasciae, plantar fasciae, trigger bands, shoulder and neck region



PERI-ACTOR SCOOP Treatment at or below the edges of bones: Front iliac crest, edge of the tibia



PERI-ACTOR® SCRAPER Treatment of all fasciae trigger band



PERI-ACTOR SPHERE Punctual treatment of fasciae: Plantar fasciae, Arch of foot, Hand, Trigger points, Cervical spine Fasciae

VIBRATION THERAPY



V10 10mm for superficlea trigger points



V25 25mm for muscle massage



V40 40mm for large muscles groups

TECHNICAL INFORMATION

2174: 230V AC 50-60Hz 350VA 2175: 120V AC 50-60Hz 400VA 2175: 100V AC 50-60Hz 400VA DIMENSIONS WXHXD: 410x1100x425mm WEIGHT: 30kg COMPRESSED AIR OUTPUT: 0.3-5.0 bar POWER INCREMENT SETTINGS: 0.1 bar PULSE AMPLITUDE: Approx 2mm at 3 bar, idling (without coupling) PULSE AMPLITUDE: Approx 2mm at 3 bar, idling (without coupling) PULSE FREQUENCY: FALCON® handpiece 1-21Hz for RPW, V-ACTOR® handpiece 1-50Hz for vibration MODE: Single pulse or continuous USER INTERFACE: 7 inch, 800x480 pixels color graphic capacitive touch screen CONFORMITY: This device complies with the applicable standards EN/IEC 60601-1, CAN/CSA-C22.2 No.601.1, UL Std. No 60601-1. CLASSIFICATION ACCORDING TO MDD: class IIa device 93/42 60601-1 LEGAL MANUFACTURER: DJ0 France SAS - CE 2797

Chothonoos



INTELECT® RPW 2

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	
2174KIT	Intelect® RPW 2 230V unit set including standard accessories	
2173KIT	Intelect® RPW 2 120V unit set including standard accessories	
2175KIT	Intelect® RPW 2 100V unit set including standard accessories	

STANDARD ACCESSORIES

PART NUMBER	DESCRIPTION				
18080	FALCON® handpiece, including 2x projectile; 2x guiding tube; 1 cleaning brush, 1x wrench; 1x allen key				
1x Power cord with EU plug for 2174KIT, 1x Power cord with North American plug for 2173KIT, 1x Power cord with Japanese plug for 2175KIT					
13-28660*	Quick start guide booklets (available in multiple languages)	28178	Transmitter Ro40		
		29724	Transmitter D20-S black		
22652	Transmitters tray	82-0274	Intelect RPW 2 pouch		
22651	Gel bottle support cup	02-02/4			
		15-1140	USB Key		
4248	1x Transmission gel, 250 ml				

OPTIONAL ACCESSORIES

PART NUMBER	DESCRIPTION			
18080	FALCON handpiece, including 2x projectile; 2x guiding tube; 1 cleaning brush, 1x wrench; 1x allen key			
29104	Revision kit short including 2x projectiles; 2x guiding tubes short; 2x o-ring set; 1x o-ring guide			
29801	A6 transmitter black	82-0274	Intelect RPW 2 pouch	
29802	T10 transmitter black	22651	Gel bottle support cup	
29729	C15 transmitter black	22652	Transmitter tray	
29728	DI15 transmitter black	4248	24x Transmission gel, 250 ml	
29726	F15 transmitter black	14820	Power cord EU	
29724	D20-S transmitter black	14679	Power cord North America	
29725	D35-S transmitter black	14821	Power cord JAPAN	
28736	D20-T transmitter black	14814	Power cord Swiss	
28178	Ro40 transmitter red	14819	Power cord Australia	
18638	R15 transmitter black	14815	Power cord UK	
29539	Atlas RPW transmitter	14816	Power cord India	
19365	V-ACTOR® HF handpiece set Including V-ACTOR®, V25 and V40 transmitters	14817	Power cord Israel	
29742	V-ALTOR [®] , V25 and V40 transmitters	14818	Power cord Denmark	
28740	V-ACTOR V25 transmitter			
28741	V-ACTOR V40 transmitter			
31800	Transmitter box including 9 transmitters (A6, T10, C15, D115, F15, D20-S, D35-S, D20-T, Ro40)			
28945	PERI-ACTOR® fascia transmitter set including 4 transmitters			
28946	SPINE-ACTOR [®] transmitter set including 3 transmitters			



DJO® Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM www.DJOglobal.eu

"Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Information related to various health, medical, and fitness conditions and their treatment is not meant to be a substitute for the advice provided by a physician or other medical professional."