

Triton Traction Unit



Ordering Information	
Part Number	Description
Triton Traction Unit	
4749	Triton Traction Unit
Standard Accessories	
48083	Patient Interrupt Switch
27465	Patient Data Cards (1x)
48144	User Manual on CD
Optional Accessories	
4766	EMG module includes:
- 48084	Patient Interrupt/EMG Indicator Switch
- 27321	EMG Leadwire
- 42192	DuraStick® Plus 5 cm Round Electrodes (4x)
- 42205	DuraStick® Premium 3 cm Round Electrodes (4x)
Technical Specifications	
Triton Traction Unit	
Mains Power:	100 - 240 V, 50/60 Hz
Power Consumption:	75 W
Current Consumption:	3.2 Amps Max
Weight:	14 kg
Dimensions*:	45 x 24 x 45 cm
Electrical Safety Class:	Class 1, Type B
Safety Tests:	IEC/UL/EN 60601-1, IEC/EN 60601-1-2, CAN C22.2 No. 601.1-M90 w/A2 - Meets MDD 93/42/EEC, CE 0413
EMG Module	
Electrical Safety Class:	Class 1, Type BF
* Dimensions are listed as L x W x H	



Features

- Static, intermittent, and cycling traction
- Independent timers allow easy customization of traction pulls in progressive, regressive, and traction phases
- Variable speed motor (30%, 50% and 100%)
- Pretension step eases the patient into traction to prevent muscle guarding and pretest patient's tolerance to traction force
- Clinical Protocols™ provided for both cervical and lumbar conditions
- Full color graphic library complete with visual tutorial of belting techniques, anatomical images and common pathologies
- Digital high resolution color touch screen that pivots 270°
- 80 User defined protocol slots
- Easily store treatment parameters, pain scales, and pain maps on Patient Data Cards
- Upgrade with the optional EMG module