

Tru Trac Traction Unit and Traction Tables



Features

Tru-Trac Traction Unit

- Static, intermittent, and cycling Traction
- User-defined hold, rest and treatment times
- 10 User-defined protocol slots to store customized settings
- Digital monochromatic touch screen interface that pivots 270°
- Easily store treatment parameters, pain scales, and pain maps on Patient Data Cards

TTET-200

The TTET-200 is a two-section, 71 cm wide table with an adjustable height from 48 cm - 96 cm and a working load of 181 kg. Optional retractable casters are available. Traction unit and traction stool sold separately. Available in all upholstery colors shown in chart.

TTFT-200

The TTFT-200 is a two-section, 71 cm wide table with a fixed height 81 cm from top to floor and a working load of up to 159 kg. Traction unit and traction stool sold separately. Available in all upholstery colors shown in chart.



Ordering Information

Part Number Description

Tru-Trac Traction Device

4779 Tru-Trac Traction Unit

Standard Accessories

48082 Patient Interrupt Switch

27465 Patient Data Cards (1x)

48079 User Manual on CD

Economy Hi-Lo Traction Table

6870 TTET-200

Optional Accessories

66898 Foot Switch

86016 Retractable Casters

Economy Fixed Height Traction Table

6840 TTFT-200

Technical Specifications

Tru-Trac Traction Device

Mains Power: 100 - 240 V, 50/60 Hz

Power Consumption: 75 W

Current Consumption: 3.2 Amps Max

Weight: 14 kg (30 lbs)

Dimensions*: 45 x 24 x 45 cm

Electrical Safety Class: Class 1, Type B

Safety Tests: IEC/UL/EN 60601-1, IEC/EN 60601-1-2, CAN C22.2 No. 601.1-M90 w/A2 - Meets MDD 93/42/EEC, CE 0413

Economy Hi-Lo Traction Table

Mains Power: 230 V, 50/60 Hz; 120 V, 60 Hz

Current: 230 V - 1.0 A; 120 V - 1.6 A

Power Input: 230 V - 230 W; 120 V - 200 W

Weight: 64 kg (140 lb)

Dimensions*: 196 x 71 x 48 to 96 cm

Electrical Safety Class: Class 1, Type B

* Dimensions are listed as L x W x H

For Color Options, please see page 5