

Ordering Information

Part Number Description

Intelect Neo Vacuum Module

70006 Intelect Neo Module Vacuum

Standard Accessories

13-8809 Vacuum Module User Manual

13-7652 Intelect Neo User Manual on CD

70040 Vacuum Module Electrode/Leadhose Kit
(includes 60 mm electrodes and sponges)

70041 Vacuum Plug Kit

Optional Accessories

70030 Leadhose Vacuum Ch 1/2 Kit

70031 Leadhose Vacuum Ch 3/4 Kit

70032 Leadhose Vacuum Ch 1/2 Kit XL

70033 Leadhose Vacuum Ch 3/4 Kit XL

70034 90 mm Vacuum Electrode Kit

70035 60 mm Vacuum Electrode Kit

70036 30 mm Vacuum Electrode Kit

70037 90 mm Vacuum Sponge Kit

70038 60 mm Vacuum Sponge Kit

70039 30 mm Vacuum Sponge Kit

70040 Vacuum Module Electrode/Leadhose Kit

70041 Vacuum Plug Kit

Technical Specifications

Mains Power: 100-240 VAC, 50-60 Hz

Weight: Module only: 2.3 kg

Dimensions*: Module only: 26 x 14 x 28 cm

Electrical Safety Class: Class I, Type BF

Safety Tests: IEC 60601-1, IEC 60601-1-2,
IEC 60601-2-10

* Dimensions are listed as L x W x H

Intelect® Neo Vacuum Module



Features

A vacuum electrode module designed specifically for the Intelect Neo Therapy system. Easily installed in the Neo's therapy cart, the vacuum module fully integrates with the Neo's advanced features such as the Clinical Protocol Setup (CPS) system.

- 2 channel module (4 electrode connections) that can support 4 channels with the addition of a second module in the Therapy System Cart
- Continuous or Pulsed vacuum in 10 steps
- Easily adjusted vacuum setting from 0 mbar to 600 mbar in steps of 60 mbar
- Quick release reservoir allows fast and clean removal of fluid
- Quiet and efficient operation
- Various sizes of Vacuum Electrode Cups available
- 2 year warranty (accessories 1 year)

