

Intelect Focus Shockwave

Features

LCD touch screen For easy control

Transportable Not limited to one room or facility

3 standoff variations For optimum therapeutic effectiveness

12.5 *cm penetration with pinpoint focus* Greater depth to treat more indications

Enhanced Energy: 0.01 – 0.55mj/mm2 Greater stimulation of metabolic change

Broad Frequency Range: 1 – 8Hz Customize treatment to patient needs


Unlike the Radial Pressure Wave, the Focus Shockwave's power allows it to disintegrate calcifications, stones, and growths as well as providing the standard benefits of shockwave therapy:

- Activation of the Healing Cascade
- Enhanced Blood Flow
- Tissue Regeneration

Indications

- Achillodynia
- Plantar Fasciitis
- Tibial stress syndrome
- Calcific tendinitis
- Patellar tendinitis
- Epicondylitis
- Metatarsalgia
- Trigger points: periarticular shoulder pain

Stand-offs and penetration depths of the F-SW handpiece







Technical Information

F-SW operating mode	F-SW: Single shock, continuous shock 1–8 Hz
F-SW energy selection	in steps from 0.01 to 0.55 mJ/mm2
Mains input voltage	100 – 240 VAC
Mains frequency	50 / 60 Hz
Mains fuse	T5AL/250 VAC
Power consumption	max. 450 VA
Ambient temperature during operation	10 – 30 °C
Ambient temperature during storage and transport	o – 60 °C frost free
Ambient pressure during storage and transport	500 – 1060 hPa
Ambient air pressure during operating	800 – 1060 hPa
Air humidity	5 – 95%, non-condensing
Control device weight	22.1 kg
F-SW handpiece weight	770 g
Housing dimensions	(W x H x D) 450 x 165 x 530 mm
Classification according to MDD	Class IIb device
Protection against the ingress of water	IPX1

Ordering Information

Part Number	Product Description
21090	Intelect F-SW Set ROW including:
- 21095	Intelect F-SW ROW unit
- 19000	F-SW Handpiece
- 4600	Water bag
- 4700	Silicone oil
- 4248	Gel bottle
- 13-00061	Operating manual
- 0.0032.012	Power cord EU
4650	F-SW cart (optional)

F-SW Hand Piece Accessorie

19100	Stand-off I (30 mm)
19200	Stand-off II (15 mm)
19300	Closing ring transparent for stand-off I and II

DJO Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM
www.DJOglobal.eu



