

Ordering Information

Part Number Description

Intellect Neo EMG + Stim Module

70004	Intellect Neo Module Stim Channels 1 & 2 + EMG
-------	--

Standard Accessories

13-8905	Stim 1/2 Module User Manual
13-7718	Modules User Manual on CD
42182	Dura-Stick® Plus 2 inch round (2 packs of 4 each)
70010	Stim Ch 1/2 Leadwire Kit
70014	Stim Ch 1/2 + EMG Leadwire Kit

Intellect Neo Stim Module

70000	Intellect Neo Module Stim Channels 1 & 2
-------	--

Standard Accessories

13-8905	Stim 1/2 Module User Manual
13-7718	Modules User Manual on CD
70010	Stim Ch 1/2 Leadwire Kit
42182	Dura-Stick® Plus 2 inch round (2 packs of 4 each)

Intellect Neo Channel 3/4 Stim Module

70003	Intellect Neo Channel 3/4 Stimulation Module
-------	--

Standard Accessories

13-8893	Stim 3/4 Module User Manual
13-7718	Modules User Manual on CD
42182	Dura-Stick® Plus 2 inch round (2 packs of 4 each)
70011	Stim Ch 3/4 Leadwire Kit

Optional Accessories for all Stim modules

70001	Intellect Neo Cart
70008	Intellect Neo Patient Remote
79977	HiVolt Probe
70012	Stim Ch 1/2 XL Leadwire Kit
70013	Stim Ch 3/4 Leadwire Kit XL
42192	Dura-Stick® Plus 5 cm round electrodes
10648	Nylatex Wrap
79967	6 x 8 cm carbon electrodes (4x)
79970	6 x 8 cm sponges (4x)

Intellect® Neo sEMG and Stimulation Modules

**Intellect Neo sEMG & Channel 1/2 Stimulation Module**

The sEMG module combines surface EMG biofeedback with 2 channel electrotherapy and the ability to trigger stim from sEMG outputs. The module fully integrates with the Neo's advanced features such as the Clinical Protocol Setup (CPS) system.

- Over 30 clinical waveforms
- Up to 4 independent electrotherapy channels (with optional channel 3/4 module)
- A comprehensive electrode placement guide with detailed descriptions and imagery
- Improved sEMG functionality offering both bar and linear graph displays
- Provides precise biofeedback of muscle activity with a sensitivity of less than 3 micro volts
- Specific muscle isolation for re-education, strengthening or rehabilitation

Intellect Neo Channel 3/4 Stimulation Module

The Channel 3/4 module gives a total of four independent electrotherapy channels and fully integrates with the Neo's advanced features such as the Clinical Protocol Setup (CPS) system.

Intellect Neo Channel 1/2 Stimulation Module

A 2 Channel Electrotherapy module (without sEMG) designed specifically for the Intellect Neo Therapy System.