



# Galaxy 3 Section High-Low



## Product features:

- Adjustable head section from +40 to -70 with gas spring
- Easy access foot switch, accessible from any position around the table 360°, for simple height adjustment
- Face hole on head section with removable face hole cover
- Adjustable leg section from 0 to 90 with gas spring
- Two-layer padding with non-flammable, bio-compatible and scratch resistant upholstery in many colours
- Powder coated, extremely stable scissor steel frame
- 2 side brackets for fixation of stabilization belts and elastic bands
- Traction support sleeves
- Personal Authorization System with 2 safety magnetic keys
- Hallotronic® System electric high/low height adjustment from 52 to 102 cm
- Non-slip rubber feet with 1 cm height adjustment allowing leveling of the table
- 2 types of frame:
  - 2 small wheels and 2 legs
  - 4 directional castors with central brake and 4 legs

## Table Articulations/Dimensions



## Ordering Information

### Part Number Description

#### Galaxy 3 Section High-Low Treatment Table

|        |  |
|--------|--|
| 3133-x | Galaxy 3 Section Electrical, All-round Switch, Castors |
| 3132-x | Galaxy 3 Section Electrical, All-round Switch          |
| 3131-x | Galaxy 3 Section Hydraulic, Castors                    |
| 3130-x | Galaxy 3 Section Hydraulic                             |

## Technical Specifications

### Electrical

|                          |                                |
|--------------------------|--------------------------------|
| Mains Power:             | 100 - 240 Volt 50/60 Hz        |
| Standby Consumption:     | max. 1 A                       |
| Weight:                  | Standard 87 kg/4 Castors 99 kg |
| Shipping Weight:         | Table Weight + 50 kg           |
| Shipping Dimensions*:    | 220 cm x 85 cm x 77 cm         |
| Lifting Capacity:        | 200 kg                         |
| Height Adjustment:       | 52 - 102 cm                    |
| Electrical safety class: | Class II, Type B               |

### Hydraulic Pump

|                       |                                |
|-----------------------|--------------------------------|
| Weight:               | Standard 88 kg/4 Castors 97 kg |
| Shipping Weight:      | Table Weight + 50 kg           |
| Shipping Dimensions*: | 220 cm x 85 cm x 77 cm         |
| Lifting Capacity:     | 200 kg                         |
| Height Adjustment:    | 51 - 96 cm                     |

## Upholstery Color Options

|   |   |   |
|---|---|---|
|  |  |  |
| Ecu-Crème<br>(x = 00)   | Beige<br>(x = 01)   | Blue<br>(x = 02)  |
|  |  |  |
| Navy Blue<br>(x = 03)   | Terra Cotta<br>(x = 04)   | Chocolate<br>(x = 05)   |
|  |  |  |
| Gray<br>(x = 06)  | Graphite Gray<br>(x = 07)   | Black<br>(x = 08)   |

NOTICE: Color charts are produced within the limitations of lithographic color printing and are intended as a general guide only.

For ordering purposes, please use model number, hyphen, upholstery color number as indicated in the chart above.