

Hydrocollator® Stationary Heating Units



Hydrocollator Accessories

Side Table Rack

Stainless steel work surface and wire shelves add efficiency to units. Racks available to fit most models (M-4, M-2).

Lotion Warmers

Keeps massage lotion and ultrasound gel at a comfortable temperature. Easily attaches to any treatment table or wall. Mini-Lotion Warmer hangs on the front of the Hydrocollator E-1 or E-2 heating units (bottles not included).

Wall Mounted Towel Rack

This stainless steel rack provides a convenient way to dry or store towels and terry covers (towels not included).



Ordering Information

Part Number Description

E-1 Stationary Hydrocollator

2102 120 V Unit

2107 220 V Unit

Standard Accessories

1006 Standard size HotPacs (4x)

E-2 Stationary Hydrocollator

2802 120 V Unit

2806 220 V Unit

Standard Accessories

1006 Standard size HotPacs (3x)

1004 Oversize HotPacs (2x)

1002 Neck Contour HotPac (1x)

Hydrocollator Accessories

Side Table Rack

4010 Side Table Rack (not shown)

4230 Extra Shelf for Side Table Rack

TM-1 Electric Lotion Warmer

4012 TM-1 Electric Lotion Warmer (110 V only)

Mini Non-Electric Lotion Warmer

22100 Mini Lotion Warmer

Dial Thermometer

4228 Dial Thermometer (not shown)

Wall Mounted Towel Rack

4016 Wall Mounted Towel Rack 5 x 84 cm

Technical Specifications

E-1 Stationary Hydrocollator

Mains Power: 110-240 V 50/60 Hz

Power Consumption: 1000 W

Weight: 7 kg

Dimensions*: 33 x 20 x 41 cm

Electrical Safety Class: Class 1, Type B

Safety Tests: EN 60601-1

E-2 Stationary Hydrocollator

Mains Power: 110-240 V 50/60 Hz

Power Consumption: 1000 W

Weight: 14 kg

Dimensions*: 38 x 33 x 51 cm

Electrical Safety Class: Class 1, Type B

Safety Tests: EN 60601-1

* Dimensions are listed as L x W x H