

Order for round knit made-to-measure
compression stockings / pantyhose

medi

medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-781
export@medi.de
www.medi.de/en

Fax order to medi • Telefax +49 921 912-781

Customer No.

Date

Commission / Made-to-measure No.

Quantity ☐ piece(s) ☐ pair(s) ☐ pants

Sender / Date / Stamp

Foot Requirements

- ☐ Foot length open toe (ℓA)
(heel to ball cA): _____ cm
- ☐ Foot length closed toe (ℓZ)
(heel to toe): _____ cm
- ☐ Soft toe

Colour

Topbands mediven elegance, mediven
comfort, mediven cotton*, mediven
plus and mediven forte

- ☐ Plain silicone topband rose (FBR)
- ☐ Silicone dot topband (NoB)
- ☐ AD: Narrow silicone dot topband
(sNoB) 2,5 cm
- ☐ Silicone dot topband motif (NoM)
- ☐ Silicone micro dot topband sensitive
(MBs)

* mediven cotton is only available with NoB
and MBs

Waist attachment

only plus and forte

- ☐ left ☐ right
- Heel ℓT: _____ cm
- Waist
circumference cT: _____ cm

Special requests

Circumferences left leg / cm		Circumferences right leg / cm	Length panty section cm
			front
		please specify length of panty section	back
		cT	ℓT
		cH	ℓH
			Lengths leg
		cK	
cG		cG	ℓG
cF		cF	ℓF
cE		cE	ℓE
cD		cD	ℓD
cC		cC	ℓC
cB ¹		cB ¹	ℓB ¹
cB		cB	ℓB
cY		cY	
cA		cA	
ℓA	Foot length	ℓA	
ℓZ		ℓZ	

Models	CCL	CCL	CCL	Lengths					
Please mark with a cross:	1	2	3	AD	AF	AG	AT	AT / U	AT / H
mediven elegance									
mediven comfort									
mediven cotton									
mediven plus									
mediven forte									
mediven angio									

All measurements are body measurements! The person placing the
order is liable for incorrect measurements. Made-to-measure
garments cannot be returned.
All made-to-measure orders must be placed in writing.

Intended purpose

Round-knitted medical compression stocking used for compression
of the lower extremities, mainly for the treatment of disorders of the
venous system.

medi. I feel better.

Customised Orders – Genumedi® Knee Bandages

Note: Please order two strap knee and ankle bandages with the medi measurement block.

Customer name

Customer number

Commission

Age

Height

☐ Male ☐ Female

Stamp

medi

medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 912 912-781
export@medi.de
www.medi.de

Measurements

☐ left ☐ right

	Circumference (cm)	Length (cm)
CF		
CE		ℓE – F* (min. 15 cm)

Colours / Beaded silicone top band

☐ Silver ☐ Sand ☐ Beaded silicone top band (round top, 5 cm)

*With additional beaded silicone band, the band is already included in the
length (i.e. ℓE – F = 15 cm gives a top band of 10-15 cm).

Minimum lengths

ℓE – F min. 15 cm
ℓE – D min. 7 cm

99.867 / 02.2019

0097241 / 10.2023