

Mobile RPW Shockwave



MOBILE RPW

Chattanooga's Radial Pressure Wave technology is now mobile. The Mobile RPW is compact and light enough to be easily carried while also delivering a host of features and benefits.

USER PROFILE

This Radial Pressure Wave clinical unit is ideally suited to the following applications:

- Orthopaedics
- Rheumatology
- Physical Medicine
- Physiotherapy
- Sport Clubs
- Acupuncture
- Aesthetic
- Veterinary
- Podiatrists
- Chiropractors

INDICATIONS

- Myofascial pain
- Muscle tendon tissue activation
- Tendon attachment & site illness
- Acupuncture Shockwave Therapy



The Mobile RPW D-Actor Applicator allows quick and efficient uses of applicator transmitters for desired therapeutic effects during treatment.


D-ACTOR® Applicator


Applications:

- Myofascial trigger points
- Other disorders of tendon insertions
- Activation of muscle and connective tissue
- Acupuncture shock wave therapy
- Cellulite treatment (stages I – III)
- Muscle and connective tissue tightening
- Elimination of lymphatic congestion
- Body shaping
- Regulation of muscle tone
- Medical massage



D-ACTOR® Transmitters





When the need for a powerful yet portable Radial Pressure Wave device is called upon in the clinic or the athletic sports field, the Chattanooga Mobile RPW is the ideal choice. Its lightweight and compact design, together with its quiet, built-in compressor makes it perfect for portable use. User friendliness is at the heart of the Mobile RPW's design. The simple, ergonomically laid out controls enable all of the key parameters to be easily changed, while an illuminated LED display clearly indicates pressure, frequency and number of pulses. Initiation of the pulses is via a trigger mechanism on the hand piece.

FACTS AND FIGURES

- Weight of the unit: Approx. 9 kg
- Dimensions (L x W x H): 32 x 32 x 21 cm
- Built in compressor
- Single and Continuous shocks
 - Frequency D-actor: 1 - 15 Hz
 - Frequency V-actor: 1 - 21 Hz
- Application pressure: 1-4 Bar, Infinitely variable with precise adjustments of 0.1 increments.
- Oscillating D-Actor
- V-Actor Hand piece option
- USB connectivity



USER INTERFACE FEATURES

- | | |
|---------------------------|------------------------|
| 1 Pressure screen 1-4 Bar | 4 Intensity adjustment |
| 2 Pulse counter screen | 5 Pulse counter reset |
| 3 Frequency screen | 6 Frequency adjustment |

ORDERING INFORMATION:

2805 Mobile RPW unit

Standard Accessories:

28699 - D-ACTOR® Hand piece applicator

17638 - R15 15mm ESWT transmitter

28724 - D20-S D-ACTOR® 20mm transmitter

28701 - Accessory kit that includes the following:

- 28737 Sealing set
- 28738 Guide tube (x2)
- 28736 Projectile (x2)

13-5182 - Conductor™ transmission gel 250 ml (8.5 oz) bottle

0938 - Clinical Guide

13-5184 - Paper User Manual

13-5183 - CD User Manual

Optional Accessories:

28745 - Carrying Bag

28725 - D35 D-ACTOR® 35 mm transmitter

28726 - F15 Focus-Lens 15 mm transmitter

28728 - DI15 Deep Impact® 15 mm transmitter

28729 - C15 CERAmA-x® ESWT 15 mm transmitter

28730 - V-ACTOR Hand Piece

DJO International

1a Guildford Business Park

Guildford

Surrey

GU2 8XG

ENGLAND

Tel: +44 (0)1483 459 659

Fax: +44 (0)1483 459 470

www.DJOglobal.eu



Together in Motion™