

| Ordering Information |  |
|----------------------|--|
| Part Number          | Description                                |
| Montane Bobath       |  |
| 3580                 | Taurus Bobath 120 V – 1 Section – 4 Wheels |
| 3581                 | Taurus Bobath 120 V – 2 Section – 4 Wheels |
| 3582                 | Taurus Bobath 230 V – 1 Section – 4 Wheels |
| 3583                 | Taurus Bobath 230 V – 2 Section – 4 Wheels |

| Technical Specifications             |  |
|--------------------------------------|--|
| Mains Power:                         | 230 V / 50 Hz, 120 V / 60 Hz   |
| Standby Consumption:                 | 3.5 W  |
| Input:                               | 230 V / 1.0 A / 50 Hz,<br>120 V / 1.9 A / 60 Hz  |
| Weight:                              | 145 kg (310lbs) plus 50 kg<br>(110 lbs) pallet for shipping                                |
| Dimensions*:                         | 201 x 120 x 49 – 101 cm  |
| Shipping Dimensions*:                | 219 x 134 x 75 cm  |
| Lifting Capacity:                    | 150 kg (330 lbs)   |
| Electrical Safety Class:             | Class II, Type B   |
| Safety Tests:                        | Meets Directive 93 / 42 / EEC,<br>Certified to CAN/CSA Standard<br>C22.2 No. 601.1-M90x/A2 |
| * Dimensions are listed as L x W x H |  |

| Upholstery Color Options  |   |   |
|---|---|---|
|  |  |  |
| Ecru-Crème<br>(x = 00)  | Beige<br>(x = 01)   | Blue<br>(x = 02)  |
|  |  |  |
| Navy Blue<br>(x = 03)   | Terra Cotta<br>(x = 04)   | Chocolate<br>(x = 05)   |
|  |  |  |
| Gray<br>(x = 06)  | Graphite Gray<br>(x = 07)   | Black<br>(x = 08)   |

NOTICE: Color charts are produced within the limitations of lithographic color printing and are intended as a general guide only.  
For ordering purposes, please use model number, hyphen, upholstery color number as indicated in the chart above.  
Example: Montane Andes 7 Postural Drainage in Chocolate - 3578-05

# Montane Taurus Bobath Tables



- Features**
- 1 and 2 section Bobath treatment table
  - Extra wide patient surface and undercarriage
  - Electric high/low height adjustment from 49 to 101 cm for easy and safe patient transfer and adjustment to the ideal working height for treatment
  - Easy access foot switch, accessible from any position around the table 360°, for simple height adjustment
  - Features smooth vertical lift with a capacity up to 150 kg
  - Available in 4 casters or 2 wheel versions
  - 2 section head support adjusts to 85° positive inclination for treatment
  - Gas springs assist in the adjustment of head section

## Table Articulations

