



Montane Colombia Fixed Tables



Features

- 2 section fixed table for examination and treatment
- Head section manually adjusts up to a 45° positive inclination
- Solid high-quality fixed table constructed to ensure maximum stability
- Paper roll holder (optional equipment)

Table Articulations



Optional paper roll holder shown

Ordering Information

Part Number Description

Colombia Fixed Tables

3595 Colombia 2 Section Fixed

Technical Specifications

Head Angle Adjustment: 0-45°

Maximum Weight Capacity: 150 kg (330 lbs)

Weight: 31 kg (68 lbs)

Dimensions*: 189 x 64 cm

Shipping Dimensions*: 190 x 65 x 15 cm

Capacity: 150 kg (330 lbs)

Paper Roll Holder for use with paper of a maximum width of 60 cm

* Dimensions are listed as L x W x H

Upholstery Color Options



Beige
(x = 01)



Black
(x = 08)