

# Montane Dolomites



## 2 Section Examination Table



### Features

- 2 section treatment table
- Includes side rails for patient safety
- Electric high/low height adjustment from 49 to 101 cm for easy and safe patient transfer and adjustment to the ideal working height for treatment
- Easy access foot switch, accessible from any position around the table 360°, for simple height adjustment
- Features smooth vertical lift with a capacity up to 150 kg
- Includes simple to use paper roll holder
- Available in 4 casters or 2 wheel versions
- Adjustable head support - positive/negative inclinations
- Face hole on head section with removable cushion
- Side brackets for stabilizing belts are standard equipment
- Gas springs assist in the adjustment of head section

### Table Articulations



### Ordering Information

#### Part Number Description

#### Montane Dolomites

3590	Dolomites 2 Section 120 V Examination Table with Side Railings and Paper Roll Holder – 2 Wheels / 2 Legs
3591	Dolomites 2 Section 120 V Examination Table with Side Railings and Paper Roll Holder – 4 Wheels
3592	Dolomites 2 Section 230 V Examination Table with Side Railings and Paper Roll Holder – 2 Wheels / 2 Legs
3593	Dolomites 2 Section 230 V Examination Table with Side Railings and Paper Roll Holder – 4 Wheels

### Technical Specifications

Mains Power:	230 V / 50 Hz, 120 V / 60 Hz
Standby Consumption:	3.5 W
Input:	230 V / 1.0 A / 50 Hz, 120 V / 1.9 A / 60 Hz
Weight:	110 kg (242 lbs) plus 30 kg pallet when shipping
Dimensions*:	201 x 69 x 49 – 101 cm
Shipping Dimensions*:	219 x 84 x 75 cm
Lifting Capacity:	150 kg (330 lbs)
Electrical Safety Class:	Class II, Type B
Safety Tests:	Meets Directive 93 / 42 / EEC, Certified to CAN/CSA Standard C22.2 No. 601.1-M90x/A2

\* Dimensions are listed as L x W x H

### Upholstery Color Options

Ecru-Crème (x = 00)	Beige (x = 01)	Blue (x = 02)
Navy Blue (x = 03)	Terra Cotta (x = 04)	Chocolate (x = 05)
Gray (x = 06)	Graphite Gray (x = 07)	Black (x = 08)

NOTICE: Color charts are produced within the limitations of lithographic color printing and are intended as a general guide only.

For ordering purposes, please use model number, hyphen, upholstery color number as indicated in the chart above.

Example: Montane Andes 7 Postural Drainage in Chocolate - 3578-05