

Cefar® Peristim Pro



Cefar Peristim Pro

Electrostimulation for treating incontinence

Cefar Peristim Pro has 9 preset programs to treat stress incontinence, mixed incontinence, urge incontinence, pain related to urinary disorders (e.g. cystitis) and sacral nerve stimulation or acupuncture point stimulation in urge incontinence treatment.

Features

- Stimulation Using Vaginal and/or Anal Probe or Surface Electrodes
- Programmable: Cefar Peristim Pro allows you to add 3 personalized programs
- Backlit Screen: Cefar Peristim Pro has a larger screen and the backlight makes it much easier to read, even if the surrounding light is poor
- Other advantages: pause function, program lock, key lock, compliance, self check
- Two year warranty



Ordering Information

Part Number Description

Cefar Peristim Pro

1480-EXP	Cefar Peristim Pro, Complete kit (EN-DE-FR-IT-NL-SE-SP) - Probe Not Included
----------	--

Standard Accessories

1240	Set of Lead Wires (2 pcs)
1566	Cefar Peristim Pro User Manual (EN-DE-FR-IT-NL-SE-SP)
1559	Cefar Bag (Black), 2ch
1495	Cefar Neck Strap
1551	Batteries, Alk 1.5V AA (2 pcs)

Optional Accessories

1570	Fast Charger 1.2V
1575	Rechargeable Batteries, Alk 1.5V AA (2pcs)

For Electrodes, please see pages 17-18

Technical Specifications

Number of Channels:	2
Number of Programs:	9 preset, 3 programmable
Constant Current:	Up to a resistance of Max 1000 Ohm/channel
Maximum Amplitude:	99.5 mA
Intensity Increments:	Minimum 0.5 mA
Max Pulse Duration/Width:	300 μ seconds
Frequency:	120 Hz
Power Source:	2x1.5V AA
Dimensions*:	120 x 30 x 50 mm
Weight:	180g, including battery
Electrical Safety Class:	Class II, Type BF
Safety Tests:	IEC 60601-1, IEC 60101-1-2

* Dimensions are listed as L x W x H

Waveforms



- Symmetrical Biphasic 100% compensated